



Reference: PTXBC005847

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase TRIM56; RNF109; RING finger protein 109; tripartite motif-containing protein 56; tripartite motif containing 56		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 56
GENE ID	81844
ORIGIN SPECIES	Human
ACCESSIONS	BC005847
VECTOR	pDONR223
NCBI SYMBOL	TRIM56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]