

GINS3-GINS complex subunit 3 (Psf3 homolog) Gene

Reference: PTXBC005879

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PSF3; DNA replication complex GINS protein PSF3; GINS complex subunit 3 (Psf3 homolog); GINS complex subunit 3		

For research use only.

GENE DESCRIPTION	GINS complex subunit 3 (Psf3 homolog)
GENE ID	64785
ORIGIN SPECIES	Human
ACCESSIONS	BC005879
VECTOR	pDONR223
NCBI SYMBOL	GINS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagaggcttattccgagtgaggctgggtgcgtggggcctgaggagaacttcttctttggacgacatcctgatgtccacgagaagctgccgtgcgcacggagaccgcatgcctgccttggcgcttcttcttgagcgggagcgcaggcgccgagactgacaacgcggtc
ccacagggttccaagcttgaaactacccttggctgctggcaaaaggacttttgacaacaagcgacggatccttctgtggaactcccaagatctaccaagagggttggaggactgtgtcagtgagatcccaatgtggtggacctccaaaaatggggcccttctacgggttggctccc
agctcctgcattttgacagtcctcagaaatgcagacatttccagtcctgtcgcagactttatcgagcttttgcgcgcacatcctcagaaatgcttaacaaagacacttcagccctggtagccaggctagacgagatggagaggggcttattcaaacagggcagaaaggact
gaatgactttcagtggtgggagaaggcaggccttctcagatcacagctccaactcgttcagaattacaagaagaaaaaactcactgatatggaagactga