



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

COX11-COX11 homolog, cytochrome c oxidase assembly protein (yeast) Gene

Reference: PTXBC005895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX11, cytochrome c oxidase copper chaperone; COX11 homolog, cytochrome c oxidase assembly protein; cytochrome c oxidase assembly protein COX11, mitochondrial; COX11P; cytochrome c oxidase subunit 11		

For research use only.

GENE DESCRIPTION	COX11 homolog, cytochrome c oxidase assembly protein (yeast)
GENE ID	1353
ORIGIN SPECIES	Human
ACCESSIONS	BC005895
VECTOR	pDONR223
NCBI SYMBOL	COX11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagggcgtctggcgtcctgagtaggggtgcgtctcttctgttgcttgccgttgatcacccctgggtctcacaacaggcgtcagagagagggtgtagagccgttcttagggccagagtggagtgggacaggaggtgccgagagaggactagtagtgcttgccctggacatggaagcgtcga
gctcttgagccccgcattccagcattgcagccgccggcgccgctaagagctcgaacccttcacacgcgcgaggaggaggagcggcgccggcgagaaacagacacctctacttagtggccgctgtccctggggcatgctggggcgctctactcgtccgtaccctttatcgg
ctctattgccagactactggacttggagatcagcagttgcaggtcatgccicagacaagattgaacaatcgttgcctgttaagatcgaatcataaattagctttaatgcagatgtgcatgcaagcttccagtgaacttttagacctcagacaacagaaatatgtgtgtgccagagagac
tgcactgcggctgttacagagttaaagctctactcagacacgaattgggaattctcacataaatattgttcaattggagcgtggacagatttacaataataatcagctctgtttttgaagcagagcgtttaatccccagaggaatagatatccagctgtttttctacatgcatcgtgaattgc
tgaagatcacaagaattatgaattgcatctctacacactctttctacacgttttgaagcaaaaggcgcaacaagttccagcttcagagataaattga

More info: Please contact us for any question: contact@proteogenix.fr