

CD58-CD58 molecule Gene

Reference: PTXBC005930

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD58 molecule; CD58 antigen, (lymphocyte function-associated antigen 3); LFA-3; LFA3; ag3; lymphocyte function-associated antigen 3; surface glycoprotein LFA-3		

For research use only.

GENE DESCRIPTION	CD58 molecule
GENE ID	965
ORIGIN SPECIES	Human
ACCESSIONS	BC005930
VECTOR	pDONR223
NCBI SYMBOL	CD58
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttctggggagcgcggggcgccctgggggtcctcagcgtgtctgcctgctgactgctttggttcatcagctgttttcccaacaatatatggtgtgtgtatgggaatgtaactttccatgtaccaagcaatgtgcctttaaaagaggtcctatggaaaaacaaaggataaa
gttgcagaactggaaaattctgaattcagagctttctcatcttttaaaatagggtttatttagacactgtgtcaggtagcctcactatctacaacttaacatcatcagatgaagatgagatgaaatggaaatgcgaatattactgataccatgaagtctcttctttatgtgtgtgagtccttccatctcc
cacactaactgtgcattgactaatggaagcattgaagtcattgcatgataccagagtattacaacagccatcgaggacttataatgtactcatgggattgtcctatggagcaatgtaaacgtaactcaaccagtatatattttaagatggaaaatgatctccacaaaaatacagtgactctta
gcaatccattatttaataacaacatcatcaatcatttgacaacctgtatcccaagcagcggtcattcaagacacagatgacattataccataccattagcagtaattacaacatgtattgtgtgtatgaatggtatgtatgcttttaa