



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

VT A1-Vps20-associated 1 homolog (S. cerevisiae) Gene

Reference: PTXBC005937

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	vacuolar protein sorting-associated protein VTA1 homolog; C6orf55; DRG-1; DRG1; HSPC228; LIP5; My012; SBP1; LYST-interacting protein 5; SKD1-binding protein 1; Vps20-associated 1 homolog; dopamine-responsive gene 1 protein; homolog of mouse SKD1-binding p		

For research use only.

GENE DESCRIPTION	Vps20-associated 1 homolog (S. cerevisiae)
GENE ID	51534
ORIGIN SPECIES	Human
ACCESSIONS	BC005937
VECTOR	pDONR223
NCBI SYMBOL	VT A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcgcttgcaaccgtgccccgtccccgcagctgcaagacatcacatcatctgaggcagcgtcaggagcatgacaaggagaccctgtgttgcttattactgtctttatgcgaatgcagactggaatgaagatcgtatgataaaacctgaatgtgcgaattttatca
aagttaatgcatcagttagaagctctaaagaagcagttgggtgataatgaagcttattactcaagaataatgtggcgtgtgccatttggagaattatcgtttgaaatgttttgtatgcagacaatgaagatcgtgtgcgacgattcacaaaacatgatcaagctcttcttactgcgaagtctttga
tagatgtctaaacagatatttggagaactcacgtgatgaagaatgtgaaccagacagatgtgcagatgaaggccaacatcacatcataattgtttaaagaatggggagactcctcagacgcctgtgtgaaftgaagaagaattgatattgaagaagaatgaagatctggagacagactctt
gcccatcagccaactcagcatcatcattctcaacttatgccccaaggcaatgcagcgaactactatgggaatcagatftccgttgcgcacactgccagctataatcacccagcagaagtgtctccagacagcaggtgtgtagcagaataatctcaaacctacacacgactatcc
tgcattatgcccccgaattcttaacataaattcccaaggcggatgtgtcttaaccaccaagaagatttctcgaagctcagaagaactcacaatatgtctgcagctgtgttcagatgaatgaagatgaacacgtctctcagacatctacaaaaggctctcagactgaatctgacacagcagaagaatga

More info: Please contact us for any question: contact@proteogenix.fr