

GNG12-guanine nucleotide binding protein (G protein), gamma 12 Gene

Reference: PTXBC005940

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-12; G-protein gamma-12 subunit; guanine nucleotide binding protein (G protein), gamma 12; G protein subunit gamma 12		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma 12
GENE ID	55970
ORIGIN SPECIES	Human
ACCESSIONS	BC005940
VECTOR	pDONR223
NCBI SYMBOL	GNG12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagcaaaacagcaagcaccaacaatatagcccaggcaaggagaaactgtgcagcagttagattagaagcctccattgaaggaataaaggtttgaaggcatcagcggacctcatgtcctactgtgaggaaatgccaggagtgacccttgcctgataggaataccaactcagaaa
accctttcaaggataaaaaaacttgcatcatcttatag