

BANF1-barrier to autointegration factor 1 Gene

Reference: PTXBC005942

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BAF; BCRP1; D14S1460; NGPS; barrier-to-autointegration factor; breakpoint cluster region protein 1; barrier to autointegration factor 1		

For research use only.

GENE DESCRIPTION	barrier to autointegration factor 1
GENE ID	8815
ORIGIN SPECIES	Human
ACCESSIONS	BC005942
VECTOR	pDONR223
NCBI SYMBOL	BANF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaacctccaaaagcaccgagactcgtggcagagcccatgggggagaaagccagtgaggcctggctgggattggtgaagtcctgggcaagaagctggaggaaaggggtttgacaaggcctatgttgccttgccagttctggtgctaaagaaagatgaagacctctccg
ggaaatggctgaagacactgtggcgccaacccaagcagtcgggactgcttcggatgccttcgagagtggtgcgacgccttctgtga