

RAG1AP1-recombination activating gene 1 activating protein 1 Gene

Reference: PTXBC005943

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAG1AP1; HsSWEET1; SCP; SWEET1; slv; sugar transporter SWEET1; RAG1-activating protein 1; RP11-540D14.5; RZPD0834D038D; recombination activating gene 1 activating protein 1; solute carrier family 50 (sugar efflux transporter), member 1; solute carrier fam		

For research use only.

GENE DESCRIPTION	recombination activating gene 1 activating protein 1
GENE ID	55974
ORIGIN SPECIES	Human
ACCESSIONS	BC005943
VECTOR	pDONR223
NCBI SYMBOL	RAG1AP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcggcgcgctttctggaactgcctcattacggagcatgcgtggtctcacccttggcatgtctccgcccgcctctcggacctcagggcacatgcgaatgacccggagtggtgacaacgtccagttctgcctttctaccacggagtcacaacacctgggctggctgagttatggg
gccttgaaggagacgggacccatcgtcgtcaacacagtggtgctgcctcagacctgtatatcttggcatatctgcattactgcctcgggaagcgtgttgcctctacagactgcaacctgctagggtccttctcctgggtatggctacttttggctcctggtacccaacctgag
gcccggcttcagcagttggcctctctgcagtgcttcaccatcagcatgtacctcaccactggctgactgctaaagtgattcaaaactaaatcaaccaatgtctcctaccactcaccattgctaccccttctcactctgcctcctgggtcctctatgggtttcgaactcagatccctat
atcatggtgtccaacttccaggaatcgtaccagctttatccgcttctggtttctggaagtaccccaggagcaagacaggaactactggctcctgcaaacctga