

UTP23-UTP23, small subunit (SSU) processome component, homolog (yeast) Gene

Reference: PTXBC005955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UTP23, small subunit processome component; UTP23, small subunit (SSU) processome component, homolog; rRNA-processing protein UTP23 homolog; C8orf53		

For research use only.

GENE DESCRIPTION	UTP23, small subunit (SSU) processome component, homolog (yeast)
GENE ID	84294
ORIGIN SPECIES	Human
ACCESSIONS	BC005955
VECTOR	pDONR223
NCBI SYMBOL	UTP23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatcacaagcagaacaatgccaagaagcatcttggtcttccgcaacaactcggagtcgcgagccgtaccagatcctgtgagcggcaccttctgtcaggcggcgctgcggggccgcatcagctgcgggagcagctgccccgctacatcaggggagacgcagctgt
gcaccacaagatgtgtgttaaaagagctagaacaattgggaaggacttatatggggcaaaactgattgcacaaaaatgccaagttcgaaattgtctcatttcaagaatgcagtgagtgatcagaatgtctgtcttccatggttgaagagggaatcctcatcattatttgggcaacacagg
atcagaattgtctgtgaaagtaaaaaagagccctggagttcctcctcatgtttattatcagaacactatggttttgacaaaccttctccaaacaattgcctttgtaaaagcagtgagtcaggtcagcttctcagtgcatgagaagaaagtatcaaacatctcaagagggaacagggttta
gtgaaaaacactgaacagagtagaagaaaaagcgcaagaaaaataagtggtcccaatcctctttagtttgaaagaaaaagaaagcaccggacacacaatcatctgtcttctgaaaagagaagaaaaagaaagaattcggacagatctaaccctaaagtactttctgagaagcag
aatgcagaaggagaatga