

CISD1-CDGSH iron sulfur domain 1 Gene

Reference: PTXBC005962

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf70; MDS029; ZCD1; mitoNEET; CDGSH iron-sulfur domain-containing protein 1; zinc finger CDGSH-type domain 1; CDGSH iron sulfur domain 1		

For research use only.

GENE DESCRIPTION	CDGSH iron sulfur domain 1
GENE ID	55847
ORIGIN SPECIES	Human
ACCESSIONS	BC005962
VECTOR	pDONR223
NCBI SYMBOL	CISD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtctgactccagttccagctacgagttgaatgatcgagcagttaccattgctgctgggacagctgcaattggttatctagcttacaagaattttatgttaaagatcatcgaataaagctatgataaaccttcacatccagaagacaaccccaagatgacatgctttgacatgga
ggatttgggagataaagctgtgtactgccgttggagggtccaaaagtccattctgtgatgggctcacacaaaacataacgaagagactggagacaatgtggccctctgatcatcaagaaaaagaaacttaa