



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF7-ring finger protein 7 Gene

Reference: PTXBC005966

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CKBBP1; ROC2; SAG; RING-box protein 2; CKII beta-binding protein 1; rbx2; regulator of cullins 2; sensitive to apoptosis gene protein; sensitive to apoptosis, zinc RING finger protein SAG, regulator of cullins 2; zinc RING finger protein SAG; ring finger

For research use only.

GENE DESCRIPTION	ring finger protein 7
GENE ID	9616
ORIGIN SPECIES	Human
ACCESSIONS	BC005966
VECTOR	pDONR223
NCBI SYMBOL	RNF7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgacgtggaagacggagaggaaacctgcgccctggcctctcactcgggagctcagggtccaagtcgggagggcgacaagatgtctccctcaagaagtggaaacgggtcgccatgtggagctgggacgtggagtgcgatactgcgccatctgcagggtgccaggtgatgga
tgcctctcttgaatctcaagctgaaacaacaagaagaaactgtttttgtctctggggagaaatataa|cattccttcacaactgtctgcattctccctgtgtgtgaaacaagaacaatcgtctccctctctcccaacaaagactgtgtgtgtcccaaaagaatgcgaagaatga