

UBE2B-ubiquitin-conjugating enzyme E2B (RAD6 homolog) Gene

Reference: PTXBC005979

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E2-17kDa; HHR6B; HR6B; RAD6B; UBC2; ubiquitin-conjugating enzyme E2 B; E2 protein; E2 ubiquitin-conjugating enzyme B; RAD6 homolog B; ubiquitin carrier protein B; ubiquitin conjugating enzyme E2B; ubiquitin-conjugating enzyme E2-17 kDa; ubiquitin-conjugat		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2B (RAD6 homolog)
GENE ID	7320
ORIGIN SPECIES	Human
ACCESSIONS	BC005979
VECTOR	pDONR223
NCBI SYMBOL	UBE2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaccccgcccgaggaggtcatgcgggattcaagcgggtacaagaggaccacctgtgggtgtcagtgccgaccatctgaaaacaacatgcagtggaatgcagttatattggaccagaaggacacctttgaagatgggtacttttaactagtaatagaattttctgaa
gaatatccaaataaacaccaactgttaggtttttatccaaaatgtttcatccaaatgtgtatgctgatggtagcatatgttttagatatccttcagaatcgatggagtcacacatatgatgtatcttctatcttaacatcaattcagtcctgctgatgaaccgaatcctaacagtcagccaatagcca
ggcagcacagctttatcaggaacaaacagagaatagagaaaagagtttcggccattgttgaacaaagctggaatgattcataa