

TCEAL7-transcription elongation factor A (SII)-like 7 Gene

Reference: PTXBC005988

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPMGp800C04260Q003; WEX5; transcription elongation factor A protein-like 7; TCEA-like protein 7; transcription elongation factor A (SII)-like 7; transcription elongation factor S-II protein-like 7; transcription elongation factor A like 7		

For research use only.

GENE DESCRIPTION	transcription elongation factor A (SII)-like 7
GENE ID	56849
ORIGIN SPECIES	Human
ACCESSIONS	BC005988
VECTOR	pDONR223
NCBI SYMBOL	TCEAL7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaaacccctgcaaaagaaacgaaggaaagccaaagtcagcgtgccaagaggaggaaaaacgccgtatggagaatttgaacgccagcaaacagaagggaatttttagacagaggctgcttcagtcctcgaagaatttaagaggacatagactataggcattttaagatga
agaaatgacaaggaggaggagatgagatggaaggtgtttggaagagataaggggtctgagaagaatttagggctctgcattctaaccataggcattctcgggaccgtccttatccatttaa