

GINS4-GINS complex subunit 4 (Sld5 homolog) Gene

Reference: PTXBC005995

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNA replication complex GINS protein SLD5; SLD5 homolog; GINS complex subunit 4		

For research use only.

GENE DESCRIPTION	GINS complex subunit 4 (Sld5 homolog)
GENE ID	84296
ORIGIN SPECIES	Human
ACCESSIONS	BC005995
VECTOR	pDONR223
NCBI SYMBOL	GINS4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgaagaagtggatttctggacaggactctgatggggtagtgaggaagtggcttaactcctgcagagctcattgaaagattggagcaggcctggaatgaaagtgttggcctgagctgctggagagcaagcctgagattgtagaatgtgtcatggaacagctggagcac
atggaagaaatctcaggagagccaaaggaggacgtgaaggtcagcatccacaaatggagatggagaggtatccgctacgtcctcagcagctacttgcgggtcgcctcatgaagatagagaagtgtttcctcatgtccttgagaaggaaaaaacacgtcctgaggggagccttcc
agcctctcgcgggaagagttggcctttgccagagagttcatggcgaacacagagtcctatctgaaaaatgtcgccttgaagcacatgcccccctaactacagaagggtgacctctttcgggcagttcccaaccagatctagatcttactgtttctgagagtgagagaacgacaagaaaac
atactggtagaaccagacacagatgagcagagggactacgtgattgacctggagaagggtcacagcacttgatccgatacaaaaccattgcacctctgttgcacatctgagctgtccagtaatttaa