

C1orf97-chromosome 1 open reading frame 97 Gene

Reference: PTXBC005997

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf97; RP11-318L16.3; long intergenic non-protein coding RNA 467		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 97
GENE ID	84791
ORIGIN SPECIES	Human
ACCESSIONS	BC005997
VECTOR	pDONR223
NCBI SYMBOL	C1orf97
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggataagaagtcactcacagaatcctgaagatgccagggctggcaaatatgaaggtaaacacaaacgaagaaagaaagcaaaacaaaccagcaccgatcccgacatagatcagtgacgtctttttcctcagatgatcctatgttctcttctcctcatcatcgtcttcaggaa
gccagacagattcaagtattgaagatgctgccaagggaataaagaagaagagaagagaaacaaataaatgggaaaaaagaaggacaaaatatga