

LLPH-LLP homolog, long-term synaptic facilitation (Aplysia) Gene

Reference: PTXBC006002

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf31; hLLP; protein LLP homolog; LLP homolog, long-term synaptic facilitation (Aplysia); human LAPS18-like protein; LLP homolog, long-term synaptic facilitation		

For research use only.

GENE DESCRIPTION	LLP homolog, long-term synaptic facilitation (Aplysia)
GENE ID	84298
ORIGIN SPECIES	Human
ACCESSIONS	BC006002
VECTOR	pDONR223
NCBI SYMBOL	LLPH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctaaaagcttacggagtaagtggaaaagaagatgcgtgctgaaaagagaaaaaatgcccaaggaggccagcaggttaaaagtattctcaactagacggtgatgttttaataaagatgttcaagagatagcaactgtggtgtacccaacccaacattgccaagag
aaaatgcaatgtgaggtaaaagatgaaaagatgacatgaaaatggagactgatattaagagaaacaaaagactcttctagaccagcatggacagtacccaatfatggatgaacaaaggcaagaaaaaggctgaaggcaagcagagaaaaagaagggaagaaagcaaa
agcagtgaaagtggcaagggtttggcctggtag