

POP4-processing of precursor 4, ribonuclease P/MRP subunit (*S. cerevisiae*) Gene

Reference: PTXBC006098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POP4 homolog, ribonuclease P/MRP subunit; RPP29; ribonuclease P protein subunit p29; hPOP4		

For research use only.

GENE DESCRIPTION	processing of precursor 4, ribonuclease P/MRP subunit (<i>S. cerevisiae</i>)
GENE ID	10775
ORIGIN SPECIES	Human
ACCESSIONS	BC006098
VECTOR	pDONR223
NCBI SYMBOL	POP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagagtgatctaccatgcatgtctcagaagaggcgaatgactccgatgtccagccttcaggagcacagcgggccgaggccttcgtgaggcccttcctgaagcgcagcacgccccgatgagcccgaggcccgaggaccagctgcagcgcaaggcggtgtcctgga
gtacttcacccgccacaagcgcaaggagaagaagaagaagcctctctgcaggcaaaaggaggagctgcggctctttgacattaaaccagagcagagatacagcctttctccctctccatgaactctggaacagtagatcaggacgtgtgcagtgggctcaagc
cagacacgcagccagatgattcaggccaagctcttaaggcagatcttcacggggctatttattcagtgacaaaatccaaatgccctcttatgtgggtattacaggaaatcctctacaggaaacaaagcacattttcaaaattatcacaaagaagaccgctgaaagtatccccaaagcta
aactgcgtgtcactgtgaaaccgatggcctttatttcacatttacgggagcaaatccagcttcggtcaagtgaacgggtctgcgaagaagttcaagcgaagggaacgattgacctgta