

CHCHD6-coiled-coil-helix-coiled-coil-helix domain containing 6 Gene

Reference: PTXBC006123

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CHCM1; Mic25; PPP1R23; MICOS complex subunit MIC25; coiled coil helix cristae morphology 1; coiled-coil-helix cristae morphology protein 1; coiled-coil-helix-coiled-coil-helix domain-containing protein 6, mitochondrial; protein phosphatase 1, regulatory s		

For research use only.

GENE DESCRIPTION	coiled-coil-helix-coiled-coil-helix domain containing 6
GENE ID	84303
ORIGIN SPECIES	Human
ACCESSIONS	BC006123
VECTOR	pDONR223
NCBI SYMBOL	CHCHD6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggagcagcgagagcagcgaggccgaggggtgtccttcggagtggagcagcgaggagcggggtccgggtgtcaggggtgtccgggtgtctgaaaaacgtggatgaaccgcatgaaggagcccagctctccacccctgtccacacatcttctacctttggccttcaagatggcaac
ttgagagccctcacaagaatccacactgccaggtcggggagcagtggtggtccagcagccctcagggtatgaaggaggggtgtcaagaggtatgaacaggagcatgtctgtatccaggataagctcttccaggtggcaagaggggaagagagaggtgccaccaagcactccaagg
catccctgccacggggaaggcagcatcagccatgagagcagaagtcagtcgggtggcggagggagctggaagcagagagcagagcagctaaagacgcccgtgacaccttctacaaggagcagctgagcgtattgagaggaagaatgctgagatgtataactgtcttcagagcaa
ttccatgagcagcctcaaatgagagcacaataaagcccgaggggtggagccgtgtgtcaggggtgagggccagattctccactgtaccgagatgccccgcatgaggtgtgtgtgtctgacgtgtgtcaaggcataccagcgtgctgagcggccccaagggt
ga