

C11orf70-chromosome 11 open reading frame 70 Gene

Reference: PTXBC006128

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf70; chromosome 11 open reading frame 70		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 70
GENE ID	85016
ORIGIN SPECIES	Human
ACCESSIONS	BC006128
VECTOR	pDONR223
NCBI SYMBOL	C11orf70
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgggcagaatcaaggcgcaggcgttcggccttgaccagaccttcagtcctatcggaaggatgattcgttatggcctttttcaagacccaaatgttatcccaattgaagttacttcagattctctggacaatggatcatattaggaactgaagtgaataaaatgaagctataaatgttc
cttgccacacagcttcaatgtcatttttcatcggttatatgatgaagatattgacgagacagtgacataattgttaaatgttttagattcttttggatccatttctcatttctgatgagttacgaagagttttgctagtggagactcagaaaaatgaaatattcagccaaccagatagagaagagttc
ctgttttgcctttcaaacatcttgccttggaggagcccttgcataatgaggatgtgattagcccatatctggaacaacaagaagcttatctataaggatctggtgagtggtcgaagaatcctcaaccaagaaaatacagattacctcttctgtctttaagtttcagcttatgattctgctggtatgtg
ctatccttcagcaagaatcatgaacagacatttctactttattgtggatcctatcaggcgtcacttcatgttttataccactgttatgtgtgggagacatgtcttaa