

TM2D3-TM2 domain containing 3 Gene

Reference: PTXBC006150

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BLP2; TM2 domain-containing protein 3; BBP-like protein 2; beta-amyloid-binding protein-like protein 2; TM2 domain containing 3		

For research use only.

GENE DESCRIPTION	TM2 domain containing 3
GENE ID	80213
ORIGIN SPECIES	Human
ACCESSIONS	BC006150
VECTOR	pDONR223
NCBI SYMBOL	TM2D3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggagggtgcgccgctgaggggctccgcgcttgtgtcgcgtgctgctctctcgcagttctgcttctgctggcggtgagcaatcgaggcgctggctcagtcataaaggatccgggccaacacgcacattcacagtagttccaggggcagcagaagtactga
aatcccacattatgtgatgaagtgtccgagcaatggtttgtgtagcaggcttcctgcagactgtatagactgcacaacaaatttctctgtacctatgggaagcctgtcacttttgactgtgcagtgaaccatctgttacctgtgtgatcaagactcaaatcccaagaacttcatttaacatg
acttgcagattttgtggcagcttctgaacagattacgaggtgtaccaactccaccagctgcatgacgggtgctctcctggcagcgctaccctgccaaactgcacgggtgcgggaccagctccactgcttgggtaacctgtactttcccaaatgctataattgcaattggactggaggtataa
gtggtctacggctctgcttaagcatcaccctggtgggtttggagcagaccgtttctacctggccagtgccgggaaggcctcggcaagctcttcagctcgggtggcctgggaatatggacgtgatagacgtcctgctcattggagttggctatgttgaccagcagatggctctttgtac
atttag