

C7orf50-chromosome 7 open reading frame 50 Gene

Reference: PTXBC006224

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C7orf50; YCR016W; chromosome 7 open reading frame 50		
For research use only.			

GENE DESCRIPTION	chromosome 7 open reading frame 50
GENE ID	84310
ORIGIN SPECIES	Human
ACCESSIONS	BC006224
VECTOR	pDONR223
NCBI SYMBOL	C7orf50
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaacagagagaaaagtctgaagtacagagaaaaagaacaaaagctgaagaaggcgtcagcagagggggccactgctggccctgaggctgcaccaagtggcggaaggagccggctccaaggcggaagctgtgctcaggccgggctggacgcagagccagagc
tgccccagaggagcagagggctcctggaagggaagctgaaaaaggaacggaagaagaggagggcagcgtctgcccgggagggcagggcctgtggccagcacccgctgccaggcgtcggggggcgaactggccctggactacctctgcagatggggccaaagcacaagaa
ctggaggttcagaagacgaggcagacgtggctcctgctgcacatgtatgacagtgcacaaggttcccgatgagcacttctccacctgctggcctacctggaggggctgcaggccggggcccgagagctgacggtgcagaaggcgaagccctgatcggggagctggatgaggagg
gctctgatccccctgcccgggaggggccagcgcacccagcagctgctgcagctgctctctag