

FAM176B-family with sequence similarity 176, member B Gene

Reference: PTXBC006241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM176B; FAM176B; C1orf78; protein eva-1 homolog B; family with sequence similarity 176, member B; eva-1 homolog B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 176, member B
GENE ID	55194
ORIGIN SPECIES	Human
ACCESSIONS	BC006241
VECTOR	pDONR223
NCBI SYMBOL	FAM176B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgccccgcgaaggacatggagttgctcagcaacagcctggctgctacggcgacatccggccaaccccgagagcttggcctctactctgtgctgggctgtctgcttggcctgtctcaccctctgctgctcgtcagcatctctgtggcgccccgccccgccccg
gggccccggtcagcgcgggacccccgagcagcaccctggagcccgaggacgacgaggacgaggagacacgggtgactcggctgggccccgacgacacgctgccgggccccgagctgtccgagagccggacggccccctcaacgtcaacgtttcagctcggcgga
ggagctggagcggcgagcggctggagggagcgcgaacggatctgctgggagatctggcgaccggcgagccggacctgctgggcacaggcacgctggggccccagccccacggccacgggcacccctggggcgcatgcaactactga