

FGF18-fibroblast growth factor 18 Gene

Reference: PTXBC006245

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FGF-18; ZFGF5; fibroblast growth factor 18		
For research use only.			

GENE DESCRIPTION	fibroblast growth factor 18
GENE ID	8817
ORIGIN SPECIES	Human
ACCESSIONS	BC006245
VECTOR	pDONR223
NCBI SYMBOL	FGF18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtattcagcgcctccgctgcactgctgtgtttacactctcgtgtgtgtcctcaggtacaggtgtgtgttgcggaggagagacgtggactccgatccacgtggagaaccagacgccccgctcgggacgatgtgagccgtaagcagctcggctgtaccagctctacagccgga
ccagtgggaaacacatccaggtcctggtggtcggcaggtacgtgcccgcggcgaggtgggacaaagtatgccagctcctagtggagacagacaccttcgtagtcaagtcggatcaaggcgaaggagacggaattctacgtgtcatgaaccgaaaggcaagctcgtggggaag
cccgatggcaccagcaaggaggtgtgtgttcacgagaaggtctggagaacaactacacggccctgatgtcggcctaagttactccggtgtgtacgtgggttcaccaagaaggcgccggaaggccccaagaccgggagaaccagcaggacgtgcatattcatgaagcgcta
cccccaaggggcagccggagcttcagaagccctcaagtacacgacgggtgaccaagaggtcccgctcgatccggccacacacccctgcctag