

COBLL1-COBL-like 1 Gene

Reference: PTXBC006264

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COBLR1; cordon-bleu protein-like 1; COBL-like 1; cordon-bleu WH2 repeat protein like 1		

For research use only.

GENE DESCRIPTION	COBL-like 1
GENE ID	22837
ORIGIN SPECIES	Human
ACCESSIONS	BC006264
VECTOR	pDONR223
NCBI SYMBOL	COBLL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacggccgaaccccgccggcagggacgccccagccaggagaaaacaaaagccaaggcaccacttctccagctgagaccaaatactgatgtctctcagctgtgattctgtagaatccactgctttcatcatggaacagaaagaaaacatgataaagacgttgaaactct
cagtggtcctacgtgggatattatcaaatctactactgttcatggcaggtacgtggagtaaatgaaataacagacccctcttattccacattgaaagtctctgtgataacgttcttgatgtttgttcttga