

APOL4-apolipoprotein L, 4 Gene

Reference: PTXBC006276

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APOL-IV; APOLIV; apolipoprotein L4; apolipoprotein L, 4; apolipoprotein L-IV		
For research use only.			

GENE DESCRIPTION	apolipoprotein L, 4
GENE ID	80832
ORIGIN SPECIES	Human
ACCESSIONS	BC006276
VECTOR	pDONR223
NCBI SYMBOL	APOL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacatcctgggtgcagctcatcacaagcgtcgggacaagcggacttttcttgggtgtgagagtgagggaagaggagctggaatgaggtgcagcaaaaccatccaggctggacagtggctggacagttccaagggctctctggggccatccctccctgtccctacagctggtaa
cagcagctccttctgtgtccattatgtgaacctgtgctgaggttctgttcttctgtgacctccagttatccacatctcagcatggccctctgtcagctggctgtctcatgactggccgaggttaagctgtgtttgtgtctga