

CHRNA1-cholinergic receptor, nicotinic, alpha 1 (muscle) Gene

Reference: PTXBC006314

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACHRA; ACHRD; CHRNA; CMS1A; CMS1B; CMS2A; FCCMS; SCCMS; acetylcholine receptor subunit alpha; acetylcholine receptor, nicotinic, alpha 1 (muscle); cholinergic receptor, nicotinic alpha 1; cholinergic receptor, nicotinic, alpha 1 (muscle); cholinergic rece		

For research use only.

GENE DESCRIPTION	cholinergic receptor, nicotinic, alpha 1 (muscle)
GENE ID	1134
ORIGIN SPECIES	Human
ACCESSIONS	BC006314
VECTOR	pDONR223
NCBI SYMBOL	CHRNA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggagccctggcctcctcctcctctctttagcctttgctcagctggcctcgtcctggcctcgaacatgagaccctgctggtggcaaaagctatftaaagactacagcagcgtggtgctggccagtggaagaccaccgccaggtcgtggaggtcaccgtggcctgcagctgatacagctcat
caatgtggatgaatgaaatcagatcgtgacaaccaatgtgctgctgaaacagcaatgggtggattacaacctaaatggatccagatgactatggcgtgtaaaaaaaatcacattcttcagaaaagatcggcggccagacctgttctatacaaatgcagatgtgacttgcattgt
caagttcaccaaaatgctcctcagtaactggccacatcacctggacacctcagccatctttaaagctactgtgagatcactcaccaccttccctttgatgaacagaactgcagcatgaagctgggcacctggacactacagcgtcgtggtgccaatcaaccgaaagcgacc
agccagacctgagcaactcatggagagcgggagtggtgatcaaggatccgggctggaagcactcctgacctattcctgctgccccgacacccctactctgacatcacctaccactctcatgcagcgcctgccccctactctcactcgtcaacgtcatccctcctgctc
ttctcctttaactggcctgtattctacctgcccacagactcaggtgggtggttgcctgactcgtcgtgta
```