

CLTB-clathrin, light chain (Lcb) Gene

Reference: PTXBC006332

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LCB; clathrin light chain B; clathrin, light chain (Lcb); clathrin, light polypeptide (Lcb)		
For research use only.			

GENE DESCRIPTION	clathrin, light chain (Lcb)
GENE ID	1212
ORIGIN SPECIES	Human
ACCESSIONS	BC006332
VECTOR	pDONR223
NCBI SYMBOL	CLTB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgatgacttggcttctctgctgctgagagcggtgccccggaggcgcgaggaggacccggcgcccttctgcccagcaggagagcgagattgcaggcatagagaacgacgagggtctcggggcacctgccggcagccatgcggccccgcgcagcgggc
cccacgagtggggtggttctgagggacatggggaccacagtcattggagatgtttcaggaggccaacggtcctgctgagtgctacgcagccattgcccaggctgacaggctgaccaggagcctgagagcatccgaagtggcgagaggagcagaggaaacggctgcaagagc
tgatgctgcatctaaggtcacggaacaggaatggcgggagaaggccaagaaggacctggaggagtggaaccagcgcagagtgacaagtagagaagaacaagatcaacaacgggcatccgaggaggcttctgtgaaggaatccaaggaggagacccagggcacagagtg
gagaaggtggcccagctatgtgacttaacccaagagcagcagcagtgcaaatgtgtcccgctgcctcgtgctcatgtccctgaagcagacgccaactgtcccgtag