



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

TRIM14-tripartite motif-containing 14 Gene

Reference: PTXBC006333

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tripartite motif protein TRIM14; tripartite motif-containing protein 14; tripartite motif containing 14		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 14
GENE ID	9830
ORIGIN SPECIES	Human
ACCESSIONS	BC006333
VECTOR	pDONR223
NCBI SYMBOL	TRIM14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr