



Reference: PTXBC006334

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERG-2; ERG2; HERG2; Kv11.2; hERG-2; potassium voltage-gated channel subfamily H member 6; eag-related gene member 2; ether-a-go-go-related protein 2; potassium channel, voltage gated eag related subfamily H, member 6; potassium voltage-gated channel, subf		

For research use only.

GENE DESCRIPTION	potassium voltage-gated channel, subfamily H (eag-related), member 6
GENE ID	81033
ORIGIN SPECIES	Human
ACCESSIONS	BC006334
VECTOR	pDONR223
NCBI SYMBOL	KCNH6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr