

VAX2-ventral anterior homeobox 2 Gene

Reference: PTXBC006336

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DRES93; ventral anterior homeobox 2		
For research use only.			

GENE DESCRIPTION	ventral anterior homeobox 2
GENE ID	25806
ORIGIN SPECIES	Human
ACCESSIONS	BC006336
VECTOR	pDONR223
NCBI SYMBOL	VAX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgatgggggcccagcgcgaccggggcccgccggcgaggctgtgtggcgccgggtggcgctgcggagacgcagcggagcggggactgcgagctgatggcggtggccacagcccaacggaggtggccgggacctcagctcagctccgcaggctc
caggagagatggagccgacagcagcggcgagccggggcccgaggcagaccactgccgcccatactgggtgcgagatgccaaagggaatttcctgcttaaggcctggacctggaccggcccaagcggacactacatcctcactgcgagcagctgta
ccgctggagatggagttccagcgtgccagtatgtgtggcgccgagcgcactgagctggcccgccagctgaacctctccgagaccagggtgaaggtctgtgttcagaacccggcaccagcagaagaagaccagcagagacctggagaagcggcgctcctcctcagcc
tccgagccttggccacctccaactctcggctgtgtggcagggccggctgtctgtgtgccagggccctagcctctgtggcgctgacccctagcctgcccaggcctacctgcccagccaggggcacctccttaggtgacccaggaactcctcccacgctcaaccgctgtcc
tcggcctcagcgtccccccactgcgccccctctgcagctgtctgtttctcggccccctctggtatctgctgcccgtacgaactgggttctcggccttcgagccatacagctggctagaacggaaagtggcgagcggccagcagctgcaagaagtaactaa