

SPEG-SPEG complex locus Gene

Reference: PTXBC006346

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPEG complex locus; APEG-1; APEG1; BPEG; CNM5; SPEGalpha; SPEGbeta; striated muscle preferentially expressed protein kinase; aortic preferentially expressed gene 1; aortic preferentially expressed protein 1; nuclear protein, marker for differentiated aort		

For research use only.

GENE DESCRIPTION	SPEG complex locus
GENE ID	10290
ORIGIN SPECIES	Human
ACCESSIONS	BC006346
VECTOR	pDONR223
NCBI SYMBOL	SPEG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagcccagctccagcagaaccgcttctctgacactggctccaaggcacccttcaaggtctcaacttatggaccagtcagtaagagaaggccaagatgtcatcatgagcatccgctgcagggggagcccaagcctgtggtctctgctgagaaaccgagcccggtgcgcccagaccagcgcgctttgcggaggaggctgaggggtggctgtgccggctgcggatcctggctgcagagcgtggcgatgctggtttctacacttgcgaagcggtcaatgagatggtgtctcgcagtgccgagcccgcttgagggtccgagggcgagtga