



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KCNG1-potassium voltage-gated channel, subfamily G, member 1 Gene

Reference: PTXBC006367

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

K13; KCNG; KV6.1; kH2; potassium voltage-gated channel subfamily G member 1; potassium channel KH2; potassium channel Kv6.1; potassium channel, voltage gated modifier subfamily G, member 1; potassium voltage-gated channel, subfamily G, member 1; voltage-g

For research use only.

GENE DESCRIPTION	potassium voltage-gated channel, subfamily G, member 1
GENE ID	3755
ORIGIN SPECIES	Human
ACCESSIONS	BC006367
VECTOR	pDONR223
NCBI SYMBOL	KCNG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr