

INVS-inversin Gene

Reference: PTXBC006370

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INV; NPH2; NPHP2; inversin; inversion of embryo turning homolog; inversion of embryonic turning; nephrocystin-2		

For research use only.

GENE DESCRIPTION	inversin
GENE ID	27130
ORIGIN SPECIES	Human
ACCESSIONS	BC006370
VECTOR	pDONR223
NCBI SYMBOL	INVS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaagtcagagaacctgtgtttgctgttcattagcatcacaagtcctgctgctgtaagagataaggggtgctctacagagctcctcgtaggaaactctgctcttaagacaaagaagatcagttgggagacaccacttatgtattgcgtgttggtgacagattggat
tgfcagatgctcttctgaagcaggagcagatgtgaataaaactgacatagccagagaacagccctccatctgcagcccagaaggcattgagaacaatcagcagcagcagatctaa