



Reference: PTXBC006377

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	c407A10.1; phosphatidylinositol N-acetylglucosaminyltransferase subunit Q; N-acetylglucosaminyl transferase component Gpi1; N-acetylglucosamyl transferase component GPI1; PIG-Q; c407A10.1 (GPI1 (N-acetylglucosaminyl transferase component)); phosphatidylin		

For research use only.

GENE DESCRIPTION	phosphatidylinositol glycan anchor biosynthesis, class Q
GENE ID	9091
ORIGIN SPECIES	Human
ACCESSIONS	BC006377
VECTOR	pDONR223
NCBI SYMBOL	PIGQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]