



Reference: PTXBC006378

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UBA7, ubiquitin-activating enzyme E1; UBA1B; UBE1L; UBE2; ubiquitin-like modifier-activating enzyme 7; UBA1, ubiquitin-activating enzyme E1 homolog B; ubiquitin-activating enzyme 7; ubiquitin-activating enzyme E1 homolog; ubiquitin-activating enzyme E1-re		

For research use only.

GENE DESCRIPTION	ubiquitin-like modifier activating enzyme 7
GENE ID	7318
ORIGIN SPECIES	Human
ACCESSIONS	BC006378
VECTOR	pDONR223
NCBI SYMBOL	UBA7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]