

TXNDC17-thioredoxin domain containing 17 Gene

Reference: PTXBC006405

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRP14; TXNL5; thioredoxin domain-containing protein 17; 14 kDa thioredoxin-related protein; protein 42-9-9; testicular tissue protein Li 214; thioredoxin (Trx)-related protein, 14 kDa; thioredoxin-like 5; thioredoxin-like protein 5; thioredoxin domain con		

For research use only.

GENE DESCRIPTION	thioredoxin domain containing 17
GENE ID	84817
ORIGIN SPECIES	Human
ACCESSIONS	BC006405
VECTOR	pDONR223
NCBI SYMBOL	TXNDC17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccgtatgaggaggtgagcgtgtccgcttcgaggagttccacggggccgtggaacagcacaatggcaagaccatttcgcctactttacgggttctaaggacgccgggggaaagctgggtccccgactgcgtgcaggctgaaccagtcgtacgagaggggctgaagcac
attagtgaaggatgtgtgttcactctactccaagtaggagaaaaagccttattggaagatccaataatgacttcagaaaaaacttgaagtaacagcagtcctacactacttaagtatggaacacctcaaaaactggtagaatctgagtgtcttcaggccaacctggtggaatgtgttctct
gaagattaa