

UBQLN4-ubiquilin 4 Gene

Reference: PTXBC006410

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A1U; A1Up; C1orf6; CIP75; UBIN; ubiquilin-4; ataxin-1 interacting ubiquitin-like protein; ataxin-1 ubiquitin-like interacting protein; ataxin-1 ubiquitin-like-interacting protein A1U; connexin43-interacting protein of 75 kDa; ubiquilin 4		

For research use only.

GENE DESCRIPTION	ubiquilin 4
GENE ID	56893
ORIGIN SPECIES	Human
ACCESSIONS	BC006410
VECTOR	pDONR223
NCBI SYMBOL	UBQLN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccagcttttgctggaagtggaaactcacaggtgcagacgccagaagtgaattcagcagcagctggagcagctcaactccatgggctcatcaatcgtgaggctaactgcaggccctgattgccagaggaggacatcaacgcagctatcgagagactctgggctccca
gctctcctaa