



Reference: PTXBC006412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MMP-15; MT2-MMP; MT2MMP; MTMMP2; SMCP-2; matrix metalloproteinase-15; MT-MMP 2; matrix metalloproteinase 15 (membrane-inserted); matrix metalloproteinase 15 (membrane-inserted); membrane-type matrix metalloproteinase 2; membrane-type-2 matrix metalloprotei		

GENE DESCRIPTION	matrix metalloproteinase 15 (membrane-inserted)
GENE ID	4324
ORIGIN SPECIES	Human
ACCESSIONS	BC006412
VECTOR	pDONR223
NCBI SYMBOL	MMP15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1