

LOC100132167-similar to hCG1993567 Gene

Reference: PTXBC006438

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to hCG1993567
GENE ID	100132167
ORIGIN SPECIES	Human
ACCESSIONS	BC006438
VECTOR	pDONR223
NCBI SYMBOL	LOC100132167
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgttactatgacactattttacagctacccaaggtgcgaaaaagaaaaaggccaagggaagaaggtggtgctaaccctgcggtcttgaaaaagcaggagaccatgaaagtggatcttccatttgagaaattggcactggacaggacatttggcattggacagacatcc
agcccaaaaggacccacttgccttgcataatggcccatctatactaggtgcagcagcagagccatcctctataagcagctag