

HDAC7-histone deacetylase 7 Gene

Reference: PTXBC006453

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	HD7; HD7A; HDAC7A; histone deacetylase 7; histone deacetylase 7A		

For research use only.

GENE DESCRIPTION	histone deacetylase 7
GENE ID	51564
ORIGIN SPECIES	Human
ACCESSIONS	BC006453
VECTOR	pDONR223
NCBI SYMBOL	HDAC7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcttctcttctcaactcagtgccatcgccctgccgagctgcaacagcagagcaaggccagcaagatcctcattgtagactgggacgtgcacatggcaacggcaccagcaaaccttctaccaagacccagtgctctacatctccctgcatgccatgacgacggcaact
tcttcccgaggagtgaggctgtgatgaggtagggtggcagcgggtgaggcttcaatgtcaatgtggcctgggctggaggtctggaccccccatggggatcctgagtacctggctgtttcaggatagtcgtgatgccatgcccgagagttctctccagacgtactctgtgtc
tgctgattgtatgctgctgaggtcaccggccccactgggtgctacatgtttctgccaatgttttgatacatgacgcagcaactgatgaacctggcaggagcgagtggtgtgctggcctggagggtggccatgacctcacgccatctgtgacgcctctgagcctgtgtgctg
ctcttctgggtaacagggtggatccctttcagaagaaggctggaaacagaacccaacctcaatgccatccgctctctggaggccgtgatccgggtgcacagtaataactgggctgcatgcagcgcctggcctctgtccagactcctgggtgcctagagtgccaggggctgacaaag
aagaagtggaggcagtgaccgcactggcgtccctctgtgggcatcctggtgctgaagataggccctggagcagctgggtgaggaggaagaacctatgaatcttaa