

C3orf26-chromosome 3 open reading frame 26 Gene

Reference: PTXBC006475

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf26; protein CMSS1; cms1 ribosomal small subunit homolog (yeast)		
For research use only.			

GENE DESCRIPTION	chromosome 3 open reading frame 26
GENE ID	84319
ORIGIN SPECIES	Human
ACCESSIONS	BC006475
VECTOR	pDONR223
NCBI SYMBOL	C3orf26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagacgatctcggagacagtggtgggagaaccagccgactggagcaggcagcagccagaagcatcagatggggaaggagaaggagacacagaagtgatgcagcaggagacagttccagttcgtaccttcagagaaaacaaacagcctaaagaatgtttttgatacaa
ccaaaggaaagaaagagaataccaccaagaccaggaaaagaagaaagaagaataactgatgttcttcgaaaaacagaacaaacagggttacctgaagacctacagaagctgatgaaggactattatagcagcagacgttggtgattgaattagaagaactgaacctgccaga
ctcctgttctcaaggccaatgattgactcacagtcttctcctacacataaaagggaattgtcctaagtggtgtaaaacttaggaagaaccacagtgagaagaatcggctcgtgatcatctgcagctcggccgtccgagccctggagctcattaggtcgtgatgacagcatcagaggag
acggcaagttataaaattatttcaaaagcacataaaaggtccaggcaggttaagttgctggagaagcgtgtgtgacacctgggttaggaactccggggagaattaaagaactgttaacaaggtggccttaatttagcccttaaaatttctgtttttgactggaactggagagatca
gaagttgaggagaatgatggacattcccagataaagaaaggaggtatcgaactctggaatgggagtgctcagctgtgcaagtcagaatcctgaaactgggccttttctaa