

C6orf168-chromosome 6 open reading frame 168 Gene

Reference: PTXBC006515

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf168; dJ273F20; failed axon connections homolog		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 168
GENE ID	84553
ORIGIN SPECIES	Human
ACCESSIONS	BC006515
VECTOR	pDONR223
NCBI SYMBOL	C6orf168
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggccaagcttccactctgacgccactgtcttggacacttggcacaggcaatgtggaccttaccaggggacaagaccgaacggctgacaaaggtgagctgatcaaccttgccatgtactgtgagaggataaggaggaaatttggccagagtggcaccacgatgatgacaata
ccatctatgagctgagagagcagcgaaggcagcaaaacccacaccccgctgctggaatttagcttttactcaaggacagagaccttgaagatgaggagcagaaaaacagttttccagaacccagacacagattttactggacactcactcttggatcggtggacatggatgact
atacagaccacgaacagtgcaagtga