



Reference: PTXBC006552

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDS028; integrin-alpha FG-GAP repeat-containing protein 2; integrin alpha FG-GAP repeat containing 2		

For research use only.

GENE DESCRIPTION	integrin alpha FG-GAP repeat containing 2
GENE ID	55846
ORIGIN SPECIES	Human
ACCESSIONS	BC006552
VECTOR	pDONR223
NCBI SYMBOL	ITFG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]