

RPL10A-ribosomal protein L10a Gene

Reference: PTXBC006791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSA19; Csa-19; L10A; NEDD6; 60S ribosomal protein L10a; NEDD-6; neural precursor cell expressed developmentally down-regulated protein 6; ribosomal protein L10a		

For research use only.

GENE DESCRIPTION	ribosomal protein L10a
GENE ID	4736
ORIGIN SPECIES	Human
ACCESSIONS	BC006791
VECTOR	pDONR223
NCBI SYMBOL	RPL10A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcagcaaaagtctctcgacaccctgtacgaggcgggtcgggaagtcctgcacgggaaccagcgcaagcgcgcaagttcctggagacgggtggagttgcagatcagcttgaagaactatgatccccagaaggacaagcgttctcgggcaccgtcaggcttaagtcactccc
cgccctaagttctctgtgtgtctctggggaccagcagcactgtgacgaggttaaggccgtggatattccccacatggacatcaggcgcctgaaaaaactcaacaagaataaaaaactgggtcaagaagctggccaagaagtatgatgcgtttttggcctcagagtctctgatcaagcaga
ttccacgaatcctcggccaggtttaataaggcaggaagtccctccctgctcacacacaacgaaacatggtggccaaagtggatgagtggaagtcacaatcaagttccaatgaagaaggtgttatgtctggtctgtagctgttggtcacgtgaagatgacagacgatgagctgtgt
ataacattcacctggctgtcaactcttgggtgctcattgctcaagaaaaactggcagaatgtccggccttatatatcaagagcaccatgggcaagccccagcgcctatattaa