

IFITM3-interferon induced transmembrane protein 3 (1-8U) Gene

Reference: PTXBC006794

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1-8U; DSPA2b; IP15; interferon-induced transmembrane protein 3; dispanin subfamily A member 2b; interferon-inducible protein 1-8U; interferon induced transmembrane protein 3		

For research use only.

GENE DESCRIPTION	interferon induced transmembrane protein 3 (1-8U)
GENE ID	10410
ORIGIN SPECIES	Human
ACCESSIONS	BC006794
VECTOR	pDONR223
NCBI SYMBOL	IFITM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatcacactgtcaaaccttctctctgtcaacagtgccagcccccacatgatgctcaaggaggagcagaggtggctgtgctggggcgccccacacccctgctccccgacgtccaccgtgatccacatccgcagcgagacccctggtcccgaccatgtgtctgg
tcctgtttcaacaccctttcatgaaccctgtgctgctgggtttatagcattcgcctactcctggaagtctaggacaggaagatggttggcgacgtgacggggcccgagcctatgctccaccgccaagtgctgaacatctgggcccctgattctgggcatcctcatgaccattctgtca
tcgtcatcccgatgctgatcttcaggcctatggatag