

AGPAT1-1-acylglycerol-3-phosphate O-acyltransferase 1 (lysophosphatidic acid acyltransferase, alpha) Gene

Reference: PTXBC006818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	1-acylglycerol-3-phosphate O-acyltransferase 1 (lysophosphatidic acid acyltransferase, alpha)
GENE ID	10554
ORIGIN SPECIES	Human
ACCESSIONS	BC006818
VECTOR	pDONR223
NCBI SYMBOL	AGPAT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccaccagccccctgcagccctgctgaccatctcaccagacacaggggaagaagcagacatcaggtgctgcactcactctgccccctggggagttggggaaagggaacgaacctggctggaggggatag