

S100P-S100 calcium binding protein P Gene

Reference: PTXBC006819

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MIG9; protein S100-P; migration-inducing gene 9 protein; protein S100-E; S100 calcium binding protein P		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein P
GENE ID	6286
ORIGIN SPECIES	Human
ACCESSIONS	BC006819
VECTOR	pDONR223
NCBI SYMBOL	S100P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggaactagacagccatgggcatgatcatagacgtctttccgataatcgggcagcagggcagcagcacccctgaccaaggggagctcaaggtgctgatggagaaggagctaccaggcttctgcagagtggaaaagacaaggatgccgtggataattgctcaag
gacctggagccaatggagatgccaggtggacttcagtgagttcatcgtgttcgtgctgcaatcacgtctgcctgtcacaagtacttgagaaggcaggactcaaatga