

CBX5-chromobox homolog 5 (HP1 alpha homolog, Drosophila) Gene

Reference: PTXBC006821

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL25; HP1; HP1A; chromobox protein homolog 5; HP1 alpha homolog; HP1-ALPHA; HP1Hs alpha; antigen p25; chromobox homolog 5 (HP1 alpha homolog, Drosophila); epididymis luminal protein 25; heterochromatin protein 1 homolog alpha; heterochromatin protein 1-a		

For research use only.

GENE DESCRIPTION	chromobox homolog 5 (HP1 alpha homolog, Drosophila)
GENE ID	23468
ORIGIN SPECIES	Human
ACCESSIONS	BC006821
VECTOR	pDONR223
NCBI SYMBOL	CBX5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaagaaacaaacggacagctgacagttctctcagagatgaggagagatgtgtgtggagaaggtgctagacagcgccgtggttaaggacaaagtggaaatctactgaagtggaaaggctttctgaggagcacaatactgggaacctgagaaaaactggattgcct
gagctaatttctgaatttatgaaaaagtataagaagatgaaggaggtgaaaaataaaccaggagagaagtcagaaagtaacaaggaggaatccaatttctcaaacagtcgcatcaaatcaaaaaaagagagagcagagcaatgatcgcctgggctttgagagagga
ctggaaaccagaaaagatcattggggcaacagattcctgtgtgatttaattgttctaatagaatggaaagacacagatgaagctgacctgttctgcaaaagaagctaattgtgaaatgtccacaattgtgatagcattttatgaagagagactgacatggcatcatatcctgaggatgcgga
aaacaaagagaaagaaacagcaagagctaa