

RBM3-RNA binding motif (RNP1, RRM) protein 3 Gene

Reference: PTXBC006825

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IS1-RNPL; RNPL; RNA-binding protein 3; RNA-binding motif protein 3; RNA binding motif (RNP1, RRM) protein 3		

For research use only.

GENE DESCRIPTION	RNA binding motif (RNP1, RRM) protein 3
GENE ID	5935
ORIGIN SPECIES	Human
ACCESSIONS	BC006825
VECTOR	pDONR223
NCBI SYMBOL	RBM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctctgaagaaggaaagctctctgtggagggtcaactttaacaccgacgagcaggcactggaagaccacttcagcagtttcggacatctctgaggtggtcgtgtcaaggaccgggagactcagcgggtccagggtttgtttcatcacttcaccaaccagagcatgttca
gttgccatgagagccatgaacggagagctctggtggtcagatccgtgtgcatcagcaggcaagtctgtcggggaaccagagaggtggtcttggggccatgggcgtggtcgagctactctagaggtggtgggaccagggtatggagtggtgaggtattatgacagtcg
acctggagggtatggatggtatggacgttccagagactataatggcagaaaccagggtgtgtatgaccgtactcaggaggaattacagagacaattatgacaactga