



Reference: PTXBC006832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SFRS6; B52; HEL-S-91; SRP55; serine/arginine-rich splicing factor 6; SR splicing factor 6; epididymis secretory protein Li 91; pre-mRNA-splicing factor SRP55; splicing factor, arginine/serine-rich, 55 kDa; serine and arginine rich splicing factor 6		

For research use only.

GENE DESCRIPTION	splicing factor, arginine/serine-rich 6
GENE ID	6431
ORIGIN SPECIES	Human
ACCESSIONS	BC006832
VECTOR	pDONR223
NCBI SYMBOL	SFRS6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]